








GIVING A MANUAL BOLUS





1. Press Bolus key. 
2. Hit OK for duration of bolus desired, up to 10 seconds. 

GIVING BOLUS WITH DOSE CALCULATION IN DRUG LIBRARY

1. Press .
2. Press .
3. Enter volume of bolus. 
NOTE: Follow same three steps, but it will ask for bolus limit vs. bolus dose.
4. Press .
5. Time to deliver bolus is displayed.
6. Press "bolus" key. .
7. Pump will give bolus and bolus will be displayed on the screen







Pump will beep at completion of bolus.

BOLUS WITH VOLUME PRE-SELECTION



1. Press .
2. Press  and set bolus limit by using  buttons.
3. Press  to confirm and start bolus.

Note: When bolus is finished a warning "beep" will sound.


CHANGE RATE, VOLUME AND TIME WHILE INFUSING



1. Press  button to access main menu.
2. Use  to select Rate, Volume, Time, etc.
3. Press left arrow button. 
4. Enter new rate/dose with  keypad.
5. Confirm with .
6. Note: Adjust rate or dose while infusing by pushing .

TITRATE WHILE INFUSING

1. Press .
2. Adjust the rate or dose, as needed.
3. Confirm with .

SHORTCUTS

Use  key to choose preferred criteria on left hand corner of display -
ie: long drug name or volume tot (shift totals).

Depending upon therapy, different options can be selected by scrolling using the  key. The choice shown on screen will remain until  is pressed for new selection.

Infusomat® SPACE

Infusion Pump System



Quick Reference



B | BRAUN
SHARING EXPERTISE

bbraunusa.com









B | BRAUN
SHARING EXPERTISE

See Operator's Manual for complete directions for use.
Rx only. ©2008 B. Braun Medical Inc. All rights reserved. IS64 12/09 JH




TO OPEN AND LATCH DOOR

1. Turn POWER  on.
2. Press  to open door.
3. After set is loaded, hold sides of pump; use both thumbs to press door firmly closed.
4. Door will engage and finish closing.


TO START INFUSION

1. Choose New Patient by pressing .
2. If prompted for Drug Library, press .
Scroll to Care area and select by pressing .
Scroll drug names by pressing  or .
Scroll drug names alphabetically by pressing  or .
3. Enter volume to be infused and rate. If prompted, enter patient weight.
4. Press  to start infusion.



TO END INFUSION

1. Press  button to stop infusion. Green LED will go off. Disconnect from patient line.
2. Press . Confirm with .

Note: Expect some resistance from the door, and continue to pull down door.

3. Open pump door. Remove tubing and close front door.
4. Press  for 3 seconds to switch pump off.











LOADING THE SET (POWER must be ON)

1. Open door .
2. Load set from right to left.
3. Press 2-hole clip over the 2 pins on the right of the pumping section.
4. Hold 2-hole clip down, align white clip to the left over ridge in center.
5. Insure the green lever is fully pushed down and the yellow light is on.
6. Press green flow clamp into green chamber on far left until green lever marked PRESS releases out. Point on clip faces towards you and points down.
Note: Yellow light turns off.
7. Thread tubing into the left and right notches to secure; close door.
8. Press  to acknowledge Original Space Line and begin.






NOTE: When trying to "power off", pump will automatically go to "Standby Mode" for 1-24 hours.

For Clinical and Technical Support, call 800-854-6851

PIGGYBACK FUNCTION

1. Stop the infusion by pressing .
2. Use up/down arrow keys  to scroll to "Special Functions"
3. Use  to enter Special Functions.
4. Use up or down arrows  to select Piggyback. Press  to select.
5. Select automatic or manual switch over to Primary. Use  to select, then press .
6. Select option to use Drug Library for Piggyback. Use  or  to answer Yes or No to Drug Library. If Yes, select Care Area and select drug.
7. Enter VTBI and rate if prompted.
8. When done, press . The word PIGGY appears in display when infusing.

RESET DATA IN STATUS MENU (Shift Totals)

1. Press  while pump is running to get Main Menu.
2. Select "Status" in Main Menu with  and press .
3. Highlight "Volume" in ml. Press  after highlighting and pump will prompt to reset data.
4. Press  to reset shift totals.