

Percussion/Vibration

Quick Reference Guide




Mode

Press **MODE** to select Percuss, Vibrate, or Both. Select **OFF** to stop therapy.



Time

Press **TIME** to choose 5, 10, or 20 minutes.



Zone

Set therapy zone by pressing **ZONE 1**, **ZONE 2**, or **ZONE 3**.

See therapy zone illustration on back of this card or hand control.



Intensity

Press **INTENSITY+** repeatedly to set therapy intensity.

1=Least Intense

8=Most Intense




Speed

Press **SPEED+** repeatedly to set therapy speed. Speed is not adjustable in Vibrate mode.

Chest Physiotherapy Explained

The purpose of chest physiotherapy (CPT) is to help patients breathe more freely and to facilitate oxygen and carbon dioxide exchange in the lungs. Percussion/Vibration therapy is used to help break up thick secretions in the lungs by repeatedly stimulating the chest wall.

Sizewise Percussion/Vibration therapy can be programmed to focus on one, two, or all three zones of the back, targeting the movement of fluid in the patient's lungs.



All therapy should be aligned with patient's lungs

This card is not meant to replace product user manual.
For complete manual, contact Sizewise.

sizewise.com 800-814-9389


Sizewise[®]
Designed to heal