

Power Drive 360™

Quick Reference Guide



This card is not meant to replace product user manual. For complete manual, contact Sizewise.


800-814-9389

sizewise.com



Sizewise®
Designed to heal

Components



- A. **Emergency Switch:** Shuts off the device in an emergency situation when it cannot be shut down in the usual manner.
- B. **Horn:** Warns others of approach.
- C. **Speed Control:** adjusts drive speed up/down.
- D. **ON/OFF Switch:** Main power switch. ON=drive wheel is lowered for transport. OFF=drive wheel is retracted.
- E. **Directional Switches:** Selects which direction the motor will drive the device.
- F. **Battery Gauge:** Indicates the remaining battery charge. Can also assist with troubleshooting of the power drive.
- G. **Control Throttle:** Activates movement.

Warning: Transport down inclines or over thresholds—including elevators—can result in damage to bed and injury to persons. Power Drive is intended for level surfaces, but can climb a 5° incline if necessary. Extra care should be taken when transitioning at top and bottom of ramps. Loss of traction may cause damage to bed and injury to persons.

Operating Instructions

1. Secure side rails to upright and locked position. Ensure power drive control is in upright position (see next page).
2. Lower bed to lowest position.
3. Unplug bed and all accessories from wall and secure cords. If installed, battery backup will allow for basic bed or blower functions for a short period of time, depending on which is plugged into battery. *Note: Bed must be unplugged prior to powering up unit for proper operation.*
4. Unlock casters at head and foot sections of bed.
5. Press power button to turn unit on.
6. Adjust speed control. Start at low speed until familiar with unit.
7. Set direction and allow for drive wheel to reposition. To change direction, come to a complete stop, press directional button, and allow drive to reposition.
8. Pull control throttle to engage drive.
9. Power unit off by pressing power button. Drive wheel will raise.
10. Lock casters, plug bed and all accessories back into wall outlet once transport is completed.

Caution: Keep hands and feet clear of wheels during bed movement. Move only at slow walking pace. Check bed and side rail width for clearance in hallway to prevent damage to bed or other objects.

Caution: DO NOT push bed with drive wheel on ground as damage to Power Drive may occur. Should you experience any type of malfunction or are unable to raise drive wheels, use emergency release button (under deck, patient right, center of bed). Push button to raise drive wheel.



To Lower/Raise Power Drive Control




1. Pull black knob



2. Lift silver handle



3. Rotate Unit to the left



4. Rest at lowest position

To return control to upright position, follow directions in reverse.

*Note: Power drive control **MUST** be in lowered position in order to remove headboard. SizeWise recommends keeping control in lowered position when not in use for quick headboard removal.*