

Operate the GLIDE™ in three steps:

- 1 Position mattress under patient
- 2 Connect air supply
- 3 Transfer patient from stretcher to bed

Note: *The GLIDE™ requires a minimum of two caregivers. Caregivers need to be carefully placed to position the patient. If soiling is possible, place protector sheet on top of mattress, dull side down, before it is positioned under the patient.*

Warning! *Never leave patient unattended when mattress is inflated and air supply is on.*

1 Position Mattress under Patient:

1. Roll mattress lengthwise towards center so that perforated side is against the bed.
2. Place mattress under patient using a “log-rolling” technique.

Note: *Ensure patient’s head is located at the same end as the “HEAD” label on mattress.*

3. Roll patient onto his or her side toward the attendant using the bed sheet to help with log roll and place rolled edge of mattress against patient.
4. Roll patient back towards opposite side to unroll the mattress (as you would when changing a bed sheet).
5. Attach the two patient-centering straps in gentle contact with the patient. Straps need not be tight. Note: Do not pull on the patient-centering straps to transfer patient.

2 Connect Air Supply:

1. Ensure the ON/OFF (I/O) switch is positioned to OFF (O).
2. Plug air supply power cord into wall electrical outlet.
3. Attach one end of the flexible hose to the air supply rubber grommet and insert other end of the flexible hose into mattress sleeve under label flap.
4. Wrap the Velcro strap around mattress sleeve.
Note: *Velcro strap must be fastened around the flexible hose, not the hose cuff.*
5. Using the air hose retention straps, secure air hose to the mattress.
Note: *The air supply unit can remain in the roller bag if needed.*

3 Transfer Patient from Stretcher to Bed:

1. Position stretcher next to patient's bed.
2. Securely engage the brakes on the Mobile Patient Support Platform.
3. Raise the stretcher side rail located opposite the patient transfer.
4. Adjust the bed or stretcher height as close to the same level as possible.

Note: *If space between patient's bed and stretcher is greater than 3", use the transfer bridge to fill gap.*

5. Before turning air supply unit on, verify the following:

- a. Side rails, accessories or sharp objects are not obstructing the path of the mattress.
- b. Air hose is free to travel with the mattress.
- c. All patient support systems are free to travel with the patient.
- d. An attendant is positioned in the direction of the patient transfer.

6. Turn ON air supply by pressing the switch to the (I) position.
7. Wait approximately 10-15 seconds for the mattress to fully inflate.
8. Once mattress is fully inflated, grasp extended pull handles of mattress while keeping your back in the neutral ergonomic upright posture.
9. With one firm continuous pull, move patient towards attendant to desired surface.
10. Turn OFF air supply by pressing the switch to the (O) position.
11. Unplug power cord from wall and air supply.

Please see operator's manual for further information at www.stryker.com.
To learn more about Agiliti solutions please go to www.agilitihealth.com.
For further questions, please call your local Agiliti office or call **1.800.847.7368**.