

# Infusomat<sup>®</sup> Space Infusion System

Safety, flexibility and convenience  
for today's busy care environments



Maximizing flexibility while safeguarding quality of care

## B. Flexible.

**B | BRAUN**  
SHARING EXPERTISE

# Infusomat® Space Infusion System

Ideal for ICU, CCU, NICU, PICU, OR, Emergency Department and Outpatient Care Settings

## Safe

- On-board drug library accommodates hard limits and soft limits; library can be customized for each care area by drug
- Clinical advisories are customizable to your needs and procedures
- Scrolling arrows eliminate decimal place and keypad bounce errors
- Alarms are audible, clearly visible and differentiated using different tones

## Flexible

- Strong, lightweight pole clamp holds up to three pumps mounted vertically or horizontally
- Pumps easily slide onto one another and lock in place
- Space Station safely and neatly houses up to four pumps

## Easy to use

- Easy-to-learn, easy-to-use intuitive user interface
- Easy navigation with large arrow buttons
- Easy to see screen display and backlit keys in any light setting
- Reduce potential line confusion with single channel design
- Easily transportable
- Interchangeable rechargeable battery



For more information, contact Customer Service at **1-800-624-2920**, or Infusion Systems Customer Support at **1-800-627-7867**.

**B | BRAUN**  
SHARING EXPERTISE

[www.bbraunusa.com](http://www.bbraunusa.com)