

THE PoseyBED 8070

DIMENSIONS:



Area: 24 square feet - Volume (Inside): 70 cubic feet
 Total Weight completely assembled: 278 lbs.
 Total Weight of the canopy only: 17 lbs.
 Electrically articulated and controlled
 Maximum Back Angle: 65° to Horizontal
 Maximum Knee Angle: 35° to Horizontal
 Maximum Knee to Leg Angle: 35° to Horizontal

BED ELECTRICAL SPECIFICATION:

- Power: 108-132 VAC
- Frequency: 60Hz
- Current Rating: 1.6 AMP (110-120V Unit)
- Classification: Class 2, Type B (Double Insulated)
- Circuit Protection: 250V, F1.6A (110-120V Unit)
- Mode of Operation: 10% Maximum duty cycle