


Bedside Guardian™ Fall Mat



- Bonded top cover helps transfer impact forces outward
- High-density composite core absorbs impact and slows forces
- Non-slip base stays where needed
- Tapered floor perimeter helps maintain staff and patient safety

The Bedside Guardian™ Fall Mat helps minimize fall-related injuries.

The OnCare[™] Bedside Guardian[™] provides a comprehensive solution that lowers patient trauma and reduces fall-related expenses. Structural features include a bonded top cover to help transfer impact forces outward, a high-density composite core to absorb impact, a non-slip base to secure the mat, and a tapered ground perimeter to help maintain staff and patient safety. An additional Nite-Glo[™] option makes the border of the mat stand out, even in a dimly lit room.



Optional Sizes

- 24" x 70"
- 36" x 70"

Visit www.oncaremedical.com to learn more.

OnCare™ Bedside Guardian™ is a proud offering of Agiliti Health, Inc. Bedside Guardian™ is a trademark of Comfortex, an Avatar Enterprises Inc. Company.

CORPORATE OFFICE

6625 West 78th Street. Suite 300. Minneapolis, MN. 55439

Toll Free: 800.847.7368 Phone: 952.893.3200 Fax: 952.893.0704 www.agilitihealth.com

©2018 Agiliti Health, Inc. BEDSIDEGUARDIAN1218

