



- Powered lift controls the elevation rate to reduce the potential for an orthostatic drop in blood pressure
- Reduces the risk of staff and patient injury by delivering the hip, knee, and ankle joint support required during the sit-to-stand motion path
- Built-in grab bars help with balance and move with the patient
- Converts facility toilets
  to taller ADA height and
  increased weight rating
  without costly construction
  and plumbing modifications
- Available in standard and bariatric model

# LiftSeat® Powered Toilet Lift™: A scientific toileting solution for safe patient handling and falls prevention

The OnCare™ LiftSeat Powered Toilet Lift features LiftTek™ core assistive technology that reduces physical force required by patients and caregivers to complete the sit-to-stand motion of getting on and off the toilet. Used as either an over the toilet or mobile bedside commode, this innovative toileting solution gradually controls the elevation of the patient, greatly reducing the risk of orthostatic hypotension that often leads to patient falls and caregiver injury. Furthermore, toileting challenges related to sternal and hip precautions are more easily followed.



"LiftSeat" has significant implications for protecting healthcare workers from muscular skeletal disorders and preventing injury to patients."

Pete Wirtz, MD Board Certified Orthopaedic Surgeon





# LiftTek™ Core Assistive Technology:

LiftTek Core Assistive Technology replicates a healthy person's motion path of rising vertically, moving forward, and tilting, all in one fluid motion, while limiting hip flexion less than 90 degrees.

"Considerably less force is required to assist patients with toileting transfer tasks using the LiftSeat\*."

John Lloyd PhD, Board Certified Ergonomist

# **Clinical Applications:**

- Bariatrics
- Orthopedic recovery and rehabilitation
- Advanced neuromuscular disease
- CVA
- Severe knee and hip arthritis
- Oncology, related to weakness
- ICU
- Geriatrics, related to weakness and balance

## Swivel locking wheels



#### Over-the-Toilet

Designed to comply with today's taller ADA toilet height and increased weight requirements, the large, swivel locking wheels of the LiftSeat lock in place over the toilet to provide a stable, height-adjustable transfer solution for the bathroom.

Accessory drop-in commode bucket



# Portable Bedside Commode

With the patent-pending drop-in commode bucket, the LiftSeat offers a portable, reliable motorized bedside commode solution.

Grab bars swing open or can be removed



#### Flat Transfer Seat

Wide, ergonomic, easy-to-clean transfer seat is built flush to the seating platform base.

### **Handheld Control**

Handheld control pendant gives caregivers complete control of the lifting process.

# **Available Configurations:**

Item Name LiftSeat 450 Model		<b>Weight Rating</b> 450 pounds	<b>Features</b> Wheels, push bar, commode bucket, AC/DC power options
LiftSeat 600 Model		650 pounds	Wheels, push bar, commode bucket, AC/DC power options
Specifications:	LS450	LS600	
Height:	35" (at push bar)	40" (at push bar)	Power Options
Seating Height:	19" to 21.5"	20" to 22"	<b>AC Power:</b> Control Box delivers 24V DC output voltage (to motor), main input voltage 110-120V AC (medically approved), connected to GFI outlet.
Width:	24"/21.5"*	32"/28"*	
Base Footprint:	26" x 27"	28" x 30"	<b>DC Power:</b> Control Box delivers 24V DC output voltage (to motor), IP 65 rated, with acoustic battery alarm @ 50% usage, high-capacity 24V DC battery, 1 x 2.9 A, IP 65 rated (medically approved).
Weight Rating:	450 pounds	650 pounds	
*Outside/inside grab bars			<b>Charging Unit:</b> Four-hour charging cycle after initial 24-hour charge, main input voltage 110-120V AC.

# Visit www.oncaremedical.com to learn more.

OnCare™ is a proud offering of Agiliti Health, Inc. LiftSeat® is a registered trademark of LiftSeat Corporation. LiftTek™ and Powered Toilet Lift™ are trademarks of LiftSeat Corporation.

#### **CORPORATE OFFICE**

6625 West 78th Street. Suite 300. Minneapolis, MN. 55439

Toll Free: 800.847.7368 Phone: 952.893.3200 Fax: 952.893.0704 www.agilitihealth.com

©2018 Agiliti Health, Inc. LIFTSEAT1218

