
Span Surfaces for Recessed Deck Beds

Specialty models to accommodate the unique challenges of recessed deck beds.




Featuring Span's proprietary Safety Edge™ two-piece bolster system, adapted to overcome the design restrictions inherent in recessed decks.



Inner bolster "cradle" deflects when patient sits on or moves toward edge of bed . . .



Compressible section to accommodate multiple lengths and foot positions



... which helps direct patient's center of gravity safely back toward center of bed.

Choose from among these leading Span surfaces







PressureGuard® Custom Care® Convertible LAL:Convertible treatment with add-on microclimate management, alternating pressure and lateral rotation (Item # CL85V29; add-on control unit item #8400 – see page 17)



PressureGuard® Custom Care® Convertible:Convertible treatment with add-on alternating pressure and lateral rotation (Item # CJ85V29; add-on control unit item #6500 – see page 16)



PressureGuard® Custom Care®: Non-powered reactive pressure redistribution support surface (Item # NP85V29 – see page 18)

