



- ▶ Available with the following therapies:
 - Low air-loss
 - Alternating pressure
 - Pulsation
 - Adjustable heel zone

- ▶ Clinical Applications:
 - Pressure ulcer injury prevention
 - Stage 1 through 4 pressure injuries
 - High risk patients requiring aggressive wound prevention
 - Excessive maceration
 - Progressive motor decline
 - Acute surgical wounds

The SkinGuard® UMS - Recessed Deck Mattress System: Easy-to-use with the features and therapies your patients need

The microprocessor controlled SkinGuard® UMS Mattress is designed to provide caregivers the therapies their patients need in a simple to use system. Innovative engineering allows the SkinGuard® UMS system to address the needs of highly complicated patients. This unique five zone mattress system is easy-to-use, ensures patient comfort, delivers reliable outcomes and is designed to operate on Recessed Deck bed frames.

SkinGuard® UMS - Recessed Deck Mattress System

Features and Specifications

Features	
True Low Air Loss	Mattress Operates On Flat & Recessed Deck Frames
Alternating Pressure	Five Zone Auto Adjust Mattress
Pulsation	Adjustable Heel Zone
Max Inflate	Air Bladder Bottom Out Protection
Lock Out	CPR Deflate
Alarms	Cable Management System
Comfort Adjustment	Breathable Non-Shear Nylon Top Sheet
Automatic Fowler Boost	Removable Air Cells Allow For Configurable Surface
High Capacity Air Output Provides 1275 LPM Air Flow	Control Unit Operates Both Standard & Bariatric Mattresses

Specifications ⁽¹⁾	Size	Type	Weight Capacity
SkinGuard® UMS (Universal Recessed Deck Mattress System)	36"W x (78" to 90")L x 8"H Adjustable Length	Recessed Deck	500 lbs.

(1) Control unit will also operate both standard and bariatric size SkinGuard® mattresses.



MADE IN USA

1395 Pico Street, Corona, CA 92881
 Customer Service 1.866.KAP.MED.1 (888.527.6331)
 Facsimile 951.340.4361
 Patented Product: US# 7454809, US# 8220090, US#7171711

Distributed By: