




umano**medical**

MEDICAL & SURGERY BED WITH SCALE AND BED EXIT DETECTION SYSTEM



10" / 25.4 cm low height with built-in bed exit detection system to minimize falls

Seamless design to maximize infection control

Easy to use and fully featured bed

LOW HEIGHT • CLEAN • EASY TO USE • FULLY FEATURED



Seamless, smooth surfaces maximizing infection control



Conveniently located USB power outlets to plug in mobile devices / laptops



Open architecture for ease of cleaning and simple installation of mattress



Clearance for lift and overbed table compatibility

MEDICAL & SURGERY BED WITH SCALE AND BED EXIT DETECTION SYSTEM

BED FEATURES

- 4 DC Motors
- Battery back-up
- 4 casters with central brake
- One hand release siderails
- Flat deck for ease of cleaning
- Open architecture
- One step Trendelenburg
- **One step vascular foot position**
- **Integrated pump holder**
- Transport handles on both head & foot boards
- 4 corner roller bumpers
- Nurse controls integrated into siderails and footboard
- Patient control integrated into the siderails
- Patient control lockout
- Boostless™ backrest system
- Backrest angle indicator
- Compliant with standard IEC 60601-2-52

OPTIONAL FEATURES

- 10" / 25.4 cm low height
- Built-in-bed scale system
- Built-in bed exit detection system
 - 2 levels of sensitivity
 - 3 levels of sensitivity
- Remote bed exit alarm connectivity
- Bed exit SideView lights
- **Screen on footboard control**
- **Trendelenburg angle displayed on footboard**
- **Auto-contour**
- **Programmable lowest height + 2" / 5 cm**
- 80" / 203.2 cm patient sleep surface
- 84" / 213.4 cm patient sleep surface
- 3 positions integrated bed extender:
 - 80" / 203.2 cm
 - 84" / 213.4 cm
 - 90" / 228.6 cm
- 5" / 12.7 cm single or dual casters (optional antistatic casters)
- 120-volt/5A insulated outlets with USB mobile device power outlets (not available for 220 V)
- **Foot end back-up control on the frame**
- CPR release handle
- Dual nightlights
- One –Touch button groups
 - 30-degree head of bed
 - One-button chair position
 - Flat bed position
 - Backlight patient position
 - Nurse call

Accessories*

- Patient helper system
- Upright oxygen bottle holder
- IV poles (Permanent & Removable)
- 6" or 10" / 15 or 25.4 cm Bed extender cushion

* To see all accessories available, visit our website

Specifications for model FL36

Overall length (with roller bumpers)	88" / 224 cm
Extended length	92" - 98" / 234 - 249 cm
Overall width	
Siderails down	39.25" / 100 cm
Siderails up	40.5" / 102.9 cm
Weight capacity	500 lb / 227 kg
Height range (to litter top)	
Low (with 5" / 12.7 cm dual casters)	10" / 25.4 cm
High (with 5" / 12.7 cm dual casters)	30" / 76.2 cm
Low (with 5" / 12.7 cm single casters)	11.5" / 29.2 cm
High (with 5" / 12.7 cm single casters)	31.5" / 80 cm
Mattress deck positioning	
Backrest	0° to 60°
Knee Gatch	0° to 30°
Trend. /Reverse Trend.	± 16°
Boostless™	4" / 10.2 cm
Patient surface	35" x 80" or 84"
	88.9 x 203.2 or 213.4 cm
Casters diameter	5" / 12.7 cm
Bed weight	375 to 438 lb / 170 to 199 kg

Electronics

Compliant with standard	IEC 60601-1 3 rd edition
Volts	120 V AC / 220-240 V AC
Ampere rating	4.25 A (9.25 A with insulated outlets) / 2.15 A
Current leakage	< 300 micro amps
Frequency	60 Hz / 50 Hz
Hospital grade plug	Yes
Battery type	Sealed rechargeable lead acid 1.2 Ah, 2 * 12 V
Scale system precision (Precise in all bed positions)	0.5% with a difference of maximum 1.1 lb / 0.5 kg

Umano Medical reserves the right to change specifications without notice.



Umano Medical

230, boul. Nilus-Leclerc, L'Islet (Québec) G0R 2C0 CANADA
 Toll free number : 1 844 409-4030
 information@umanomedical.com
 www.umanomedical.com