

8060



Posey StaySafe™ Bed 8060*

U.S. PATENTS 7,047,991 AND 7,383,849

INDICATIONS: Patients at extreme risk of a serious injury from a fall, or unassisted bed exit.

- Enclosed Patented Mattress Compartment helps reduce the risk of patient entrapment and is expandable to accommodate 10-inch (25.4 cm) air mattress types.
- Water-repellent, skin friendly, and durable nylon canopy.
- Foam padded A-Frame holds the nylon canopy in place.
- Four access panels for patient care from all four sides.
- Quick-Release Buckles lock the zipper pull-tabs for added patient security.
- Ten Access Ports: eight tube ports and two drainage ports.
- IV Post Mount Brackets firmly attach the canopy to the hospital bed.

Cat. Description

- 8060** Posey StaySafe Bed 8060
- 8091D** Posey StaySafe Bed 8060 In-Service and Home Care DVD



Expandable Patented Mattress Compartment



Quick-Release Buckles



Access Panels on all four sides



Ten Access Ports



A-Frame with Foam Canopy Pads



*Hospital bed and mattress not included