



EZ Way, Inc.

eZ Way Smart Stand™

Introducing EZ Way's new Smart Stand™ with increased intelligence and upgraded features.

DUAL PURPOSE SIT-TO-STAND IN WEIGHT CAPACITIES OF 400, 500 & 800 LBS.

EZ Way's new proprietary circuit board system includes added features to provide customers with added intelligence and safety features including:

- **Backlit LCD screen**
During operation, EZ Way's LCD screen is now illuminated for ease in viewing lift readout information.
- **Usage time captured**
Lift usage time is captured in minutes, great for tracking frequency of use.
- **Inspection reminder**
An optional feature, helpful to ensure preventative maintenance compliance.
- **Overload protection**
Should the actuator be overloaded due to lifting excess weight, lifting against an obstruction, or mechanical failure, "Actuator Overload" will appear on the LCD screen, the stand will cease lifting, and will only lower. A great feature to protect both patients and the stand.
- **Weight Overload**
Available in units with a scale, this feature will display "Weight Overload" on the LCD screen and will only allow the stand arms to lower.
- **Weight recall**
On units with scale, the last weight taken is retained until a new weight is taken, allowing caregivers easy reference in busy working environments.
- **New safety tab design**
Strong, durable, wire clip design.
- **Easily adjustable shin pad**
Spring loaded pin is easily adjusted locking shin pad into place.



400, 500 & 800 lb. Capacities



Ambulation Aid

EZ Way Smart Stand™ Applications:

- Designed to facilitate toileting, changing of briefs and conducting pivot transfers for weight bearing patients.
- Transfer patients to/from bed, commode, wheelchair or chair.
- Ambulation feature for therapy applications.
- Patients allowed to stand safely and securely, increasing muscle strength and circulation.

eZ Way Smart Stand™ features & specifications



- Removable and adjustable footplate
- Adjusts to 3 heights for various sized patients
- Footplate removes for ambulation
- Shin pad adjusts for various sized patients
- Shin pad safety strap



HARNESSES from Small to XXXL.



STAND SUPPORT STRAP for vertical standing assistance and support.



SEAT STRAP for added safety and patient piece of mind.



PATENTED WALKER HANDLES for ambulation and standing support (available in 'long' version).



LCD READOUT
Battery gauge, function indicator, dual up/down control, digital scale readout.



NEOPRENE HARNESS
Gripping feature to help harness stay in place.



WIPEABLE HARNESS
For infection control. Less frequent need to launder.



DISPOSABLE HARNESS
Single patient use for infection control.

ADDITIONAL FEATURES:

- 10 Year Frame / 3 Year Component Warranty
- Highest Weight Capacities in the Industry
- 3 Position Base Opening For Function Around Wheelchair, Chair And Commode
- Battery Powered; Always Available For Use
- Pneumatic Hand Control
- Digital Scale Available
- Low-base, 2" And 5" Casters Available For Low Clearance Of Special Beds Or Ease Of Function On Carpeted Flooring
- Redundant Controls
- Emergency Stop Safety Control
- Emergency Lowering Option
- Anti-pinch Feature
- 24-hour Service Hotline
- Includes 2 batteries, remote charger and one harness

SPECIFICATIONS	400 lb.	500 lb.	800 lb.
Battery Power:	2 - 24 Volt lead acid	2 - 24 Volt lead acid	2 - 24 Volt lead acid
Remote Charging Unit:	1 - 120V UL listed	1 - 120V UL listed	1 - 120V UL listed
Base Dimensions w/3" Standard Caster:	24" w x 39" l x 4.5" h	26" w x 39" l x 4.5" h	33" w x 52" l x 4.5" h
Base Height w/2" Caster:	3.25"	3.25"	3.25"
Harness With Unit:	Size Medium	Size Medium	Size X-Large
Harness Bag:	Standard	Standard	Standard
Digital Scale:	Optional	Optional	Optional
Stand Accessories:	Optional	Optional	Optional