

BARIATRIC COMMODES: 700-1000 LB. WEIGHT CAPACITY



Wheeled Shower Commode Chair

Model Series 5200, 5215, 5225

- **Weight Limit: 750 pounds**
- Seat Widths: 24", 26" or 28"
- Seat Depth: 20"
- Seat to Floor Height: 21"
- Arms: Removable arms
- Cushioned seat and back for added comfort and support
- Hook-on removable footrests
- Corrosion resistant
- Hygiene access
- Rear anti-tip attachments



Bed Side Commode

Model Series 5232, 5233

- **Weight Limit: 700 - 1000 pounds**
- Seat Widths: 20", 24", 26", 30", 36"
- Seat Depth: 20"
- Seat to Floor Heights: 20" or adjustable 20", 21", 22", 23", 24"
- Arms: Removable drop down arms
- Cushioned seat
- 36"-wide model rated at 1000 lbs.
- Optional wheel kit available



Commode

Model Series 5135

- **Weight Limit: 850 pounds**
- Seat Width: 26"
- Seat Depth: 19"
- Seat to Floor Height: 21"
- Arms: Dual swing away arms

WARRANTY:

Three-year warranty
on the frame.
One-year warranty
on all other parts.