



- ▶ Specially engineered foam layers provide targeted support and pressure redistribution
- ▶ Ultra-firm edges facilitate safer, more secure transfers
- ▶ Unique heel suspension system creates a low pressure environment
- ▶ Comfortable therapeutic surface that meets the needs of heavier patients both standard and bariatric

Rest-Q® GS™ Pressure Redistribution Mattress: Therapeutic care for a wide range of patients

The OnCare™ Rest-Q® GS™ combines the Rest-Q's famous body absorbing comfort and pressure redistribution with therapeutic weight limit to meet the needs of today's patient population. The ultra-soft portions of the mattress allow the patient's body areas that are more prone to breakdown to virtually float between the mattress' supportive materials. When compared to costly therapeutic surfaces, the Rest-Q GS has demonstrated superior pressure redistribution, comfort and therapeutic weight limit.*



- ▶ Optional secure edges help retain patients on the sleep surface

*Test data and studies available upon request.



Rest-Q® GS™ Mattress Features and Specifications



Standard Features

- 1 Vapor proof, shear-reducing top cover
- 2 High-resiliency top cushion
- 3 Upper torso support
- 4 Pressure redistribution centerline
- 5 Extra support beneath the hips
- 6 Firm base layer provides support and helps prevent bottoming out
- 7 Heel suspension system
- 8 Ultra-firm perimeter support



Product	Sizes*	Weight Capacity	Warranty	Optional Features		
				Secure Edge	CFR 1633 Fire Blanket	Nylon Web Handles
Rest-Q® GS™	35" W x 80" L x 7" H	500 lbs.	5 years	✓	✓	✓
	35" W x 84" L x 7" H					
	42" W x 80" L x 7" H					
	48" W x 84" L x 7" H					
Rest-Q® GS™ Plus	35" W x 80" L x 8" H	800 lbs.	1 year	N/A	✓	✓
	35" W x 84" L x 8" H					
	42" W x 80" L x 8" H					
	42" W x 84" L x 8" H					
	48" W x 80" L x 8" H					
	48" W x 84" L x 8" H					

*Customized sizes and removable bolsters available.
Note: Specifications subject to change without notice.

Visit www.oncaremedical.com to learn more.

OnCare™ is a proud offering of Agiliti Health, Inc. Rest-Q® is a registered trademark of Comfortex, an Avatar Enterprises Inc. Company

CORPORATE OFFICE

6625 West 78th Street. Suite 300. Minneapolis, MN. 55439

Toll Free: 800.847.7368 Phone: 952.893.3200 Fax: 952.893.0704 www.agilitihealth.com

©2018 Agiliti Health, Inc. RESTQI218

