

# Spirit Select® Bed

The Intelligent Hospital Bed

**With the right programs and processes in place, the Stryker Spirit Select Bed is designed to help keep both patients and staff safe.**


Spirit Select is a low specialty acute care bed with a low bed height of 10.75 inches from floor to litter deck. Low height enhances ingress and egress from the bed by adding additional stability. The Select bed can also help improve staff efficiency and patient outcomes with Nurse Call Connectivity, SD Card Update, 32-Character LED Display, Auto Functions and Watchdog® Bed Status Indicator.



**Short hand rail allows patient support during transfers.**



**Bed Exit Alarm helps protect staff and patients from common alarm hazards such as false alarms, alarm desensitization and nonrestoration of settings.**



**Any or all patient controls can be locked or removed to restrict use with the Patient Lockout buttons.**

**10.75-inch low bed height allows patient to enter and exit with stability**

**Programmable night light adds an additional element of safety during patient transfers**

**Built-in smart scale allows for accuracy +/- 3% and is equipped with a weight change alarm**

**Nurse call connectivity provides better staff and patient communication**

# Stryker Spirit Select<sup>®</sup> Bed



## Base Model

- Battery back-up
- Patient lock-outs
- Central Lock-and-Steer Brake System
- Footboard control panel
- Auto contour
- Low chair positioning
- Trendelenburg / reverse Trendelenburg
- Emergency CPR — automatic
- Emergency CPR — manual
- High quad rails
- Total bed lock-out
- Moulded plastic headboard/footboard
- Embedded patient controls in siderails
- Integrated nurse call button
- Backlit staff controls
- Power status indicator
- Level-all bed reset
- Built-in scale
- Integrated single-zone bed exit alarm
- Programmable night light
- 4-inch tool-less bed extension

## Middle Model

(includes all features listed above, plus:)

- 120 V AC accessory outlet
- Obstruction sensors

## Loaded Model

(includes all features listed above, plus:)

- Watchdog<sup>®</sup> Bed Status Indicator

## Specifications

|                              |                         |
|------------------------------|-------------------------|
| Model Number                 | 5700                    |
| Overall Length               |                         |
| Minimum                      | 91.75" (233 cm)         |
| Maximum                      | 95.75" (243 cm)         |
| Overall Width                | 40" (101.6 cm)          |
| Weight Capacity              | 500 lb. (227 kg)        |
| Height Range (to litter top) |                         |
| High                         | 34" (86 cm)             |
| Low                          | 10.75" (27 cm)          |
| Litter Positioning           |                         |
| Knee Gatch                   | 0°–25°                  |
| Trend./Reverse Trend.        | 14° (+/-)               |
| Patient Surface              |                         |
| Minimum                      | 35" x 80" (89 x 203 cm) |
| Maximum                      | 35" x 84" (89 x 213 cm) |
| Caster Diameter              | 5" (12.7 cm)            |

## Electronics

|       |          |
|-------|----------|
| Volts | 120 V AC |
|-------|----------|

*Stryker reserves the right to change specifications without notice.*