

Stryker S3™ Med/Surg Bed

Medical

Siderails are BackSmart™ designed with one-hand operation and an intermediate position for patient assistance during bed entry and exit

iBed™ Awareness and three-zone Chaperone® Bed Exit System monitor bed status and patient position to help reduce falls

One-touch in-bed scale provides accurate weight with the patient and the bed in any position

Open architecture of S3 accepts any standard support surface, any form of traction, and is compatible with the vast majority of nurse call and networking systems

Exclusive retractable frame allows proximal location to remain essentially unchanged when the head of the bed is raised, reducing the need for patient repositioning

Stryker S3 Signature Series bed can be ordered with customized colors, wood laminate inserts, and insignia, logo or award images to support the identity of your health-care organization.



Stryker S3 is Safe, Simple and Secure – a Truly BackSmart Solution

The Stryker S3 is truly safe and secure as it monitors patient position in three zones and 30-degree Head Of Bed angle to reduce the risk of ventilator-associated pneumonia. The Stryker S3 offers one-hand operation of siderails, center-mounted brake, one-touch in-bed scale and superior ergonomics to reduce bending and twisting that can injure caregivers.



Siderails have an intermediate support position to allow patients to help themselves enter or exit the bed — reducing risk of falls and injury to the caregiver.



iBed Awareness and the three-zone Chaperone Bed Exit System monitor patient position and bed status. Choose local, remote, audible and visual signals.



One-touch in-bed scale provides accurate weight with the patient and the bed in any position.

Stryker S3

Med/Surg Bed



stryker[®]

Medical

3800 E. Centre Ave.
Portage, MI 49002 U.S.A.
t: 269 329 2100 f: 269 329 2311
toll free: 800 787 9537

www.stryker.com

Standard Features

- Centrally-located brake pedal
- 4-wheel, steel brake with swivel lock system
- 6-inch (15 cm) casters
- One-touch 30 degree Head Of Bed angle and lock
- 500 lb (227 kg) weight capacity
- One-hand, intermediate siderails
- Retractable frame
- Open architecture — flat litter deck
- Fowler angle gauge
- Four IV pole and traction equipment sockets
- Four foley bag hooks
- Steer control
- Integrated pump holder
- Foot prop rod
- Downward motion interrupt
- CPR release and removable headboard
- Dual pedestal lift design
- Three-sided bed motion controls
- One-button cardiac chair position
- Foot end lock-out controls
- Foot end digital Trend. and fowler angle display
- Trend./Reverse Trend. controls
- Neutral accent panels
- Night light

EX Model: All Standard Features Plus

- Chaperone Bed Exit System
- One-touch scale with isolated foley bag hooks
- Three-sided gatch/fowler controls

ZX Model: EX Features Plus

- Zoom[®] Motorized Drive System

Optional Features

- One-touch, in-bed scale system
- Chaperone Bed Exit System
- Chaperone with Zone Control[®] (three levels of sensitivity)
- iBed Awareness smart bed monitor
- Siderail communication package with outside gatch/fowler (Includes: nurse call, speaker, room/read light, TV/radio volume control)
- Siderail communication package with outside gatch/fowler and Smart TV (Includes: nurse call, speaker, room/read light, TV/radio volume control) and Smart TV (channel up/down, mute)

Specifications

| | |
|---------------------------------------|--------------------------------------|
| • Model Number | 3002S3 |
| • Overall Length | 93" (236 cm) |
| • Overall Width | |
| Siderails Up | 42.5" (108 cm) |
| Siderails Down | 40" (102 cm) |
| • Weight Capacity | 500 lbs (227 kg) |
| • Height Range (to litter top) | |
| High | 30" (76 cm) |
| Low (standard with 6" casters) | 16" (41 cm) |
| • Litter Positioning | |
| Backrest | 0° – 60° |
| Knee Gatch | 0° – 40° |
| Trend./Reverse Trend. | +12°/-10° (+10°/-10° with ZX Option) |
| • Patient Surface | 35" x 84" (89 x 213 cm) |
| • Caster Diameter | 6" (15 cm) |

Electronics

| | |
|------------------------------|--------------------|
| • Agency Approvals | C-US UL |
| • Volts | 115 VAC |
| • Ampere Rating | 8 A |
| • Current Leakage | < 100 microamperes |
| • Frequency | 60 Hz |
| • Hospital Grade Plug | Yes |

Stryker reserves the right to change specifications without notice.

- Mattresses: management and prevention
- 110 volt outlet power supply (not available with Zoom)
- Adapter bracket for patient helper
- Bed extender with pad
- CPR Board — removable
- Defibrillator tray/charting surface
- Upright oxygen bottle holder
- Roller bumper set (head end)
- WallSaver™ Breakaway Wall Connector
- Permanent folding IV pole (head end)
- Dual head-end 2 stage IV poles

- IV pole — removable
- Pendant control
- Pillow speaker interface
- Traction sleeve set (choose flush, 4" x 1/2", 4" x 3/4", 8" x 1/2", 8" x 3/4")
- Wood laminate options
- Custom accent color options



Mixed Sources
Product group from well-managed
forests, controlled sources and
recycled wood or fibre
Cert no. SW-COC-1530
www.fsc.org
© 1996 Forest Stewardship Council

Copyright © 2008 Stryker
PC/SCG-7.5M-0508, 07-SMD-2002
Mkt Lit-224-041708 Rev A