

Mattress System



- ▶ Standard and Bariatric models available with the following therapies:
 - Low air-loss
 - Alternating pressure
 - Pulsation
 - Adjustable heel zone
- ▶ Clinical Applications:
 - Stage 1 through 4 pressure injuries
 - High-risk patients requiring aggressive wound prevention
 - Excessive maceration
 - Acute surgical wounds

The SkinGuard[®] Mattress: Easy-to-use with the therapies your patients need

The microprocessor controlled SkinGuard[®] mattresses are designed to provide caregivers the therapies their patients need in a simple to use system. Innovative engineering allows the SkinGuard[®] systems to address the needs of highly complicated patients; including standard, bariatric, length extension, and heel pressure injuries. These unique five zone mattress systems are easy-to-use, ensure patient comfort and deliver reliable clinical outcomes.

SkinGuard® Mattress System Features and Specifications

Features	
Low Air-Loss	Five Zone Auto Adjust Mattress
Alternating Pressure	Adjustable Heel Zone
Pulsation	Air Bladder Bottom Out Protection
Max Inflate	CPR Deflate
Lock Out	Cable Management System
Alarms	Breathable Non-Shear Nylon Top Sheet
Comfort Adjustment	Output Provides 1275 LPM Air Flow
Automatic Fowler Boost	Control Unit Operates Both Standard & Bariatric Mattresses

Specifications	Size	Type	Weight Capacity
SkinGuard® Standard Replacement Mattress	36"W x 80"L x 8"H	Flat Deck	500 lbs.
	36"W x 84"W x 8"H	Recessed Deck	500 lbs.
SkinGuard® Expandable (EL) Replacement Mattress	36"W x 80"/84"/88" L x 8" H	Flat Deck	500 lbs.
SkinGuard® Bariatric Expandable (XMS) Replacement Mattress	Width expands from 36" to 42" or 48" and Length expands from 80" to 84" or 88"	Flat Deck	1000 lbs.



MADE IN USA

1395 Pico Street, Corona, CA 92881
 Customer Service 1.866.KAP.MED.1 (888.527.6331)
 Facsimile 951.340.4361
 Patented Product: US# 7454809, US# 8220090, US#7171711

Distributed By: