

Mattress System



- ▶ Standard and Overlay models available with the following therapies:
 - Low air-loss
 - Alternating pressure
- ▶ Clinical Applications:
 - Pressure ulcer injury prevention
 - Overlay Mattress: Stage 1 through 2 pressure injuries
 - Mattress: Stage 1 through 4 pressure injuries
 - High risk patients requiring aggressive wound prevention

The SkinGuard® Float Mattress System: Easy-to-use with the therapies your patients need

The microprocessor controlled SkinGuard® Float is designed to provide caregivers the therapies their patients need in a simple to use system. Innovative engineering allows the SkinGuard® Float system to address the needs of complicated patients. This unique five zone mattress system is easy-to-use, ensures patient comfort and delivers reliable clinical outcomes.

SkinGuard® Float Mattress System

Features and Specifications

Features	
Low Air-Loss	Five Zone Auto Adjust Mattress
Alternating Pressure	Air Bladder Bottom Out Protection
Max Inflate	CPR Deflate
Lock Out	Cable Management System
Alarms	Breathable Non-Shear Nylon Top Sheet
Comfort Adjustment	Removable Air Cells Allow For Configurable Surface
Automatic Fowler Boost	Extension Port - Lengthen From Standard 80" to 84" to 88" ⁽²⁾
Stays Inflated Without Power	Quiet 50 liter Compressor Provides 42-52 LPM Air Flow

Specifications ⁽¹⁾	Size	Type	Weight Capacity
SkinGuard® Float Standard Replacement Mattress	36"L x 80"W x 8"H	Flat Deck & Recessed Deck	500 lbs.
	36"L x 84"W x 8"H	Flat Deck & Recessed Deck	500 lbs.
SkinGuard® Float Overlay Mattress	36"L x 80"W x 3"H	Flat Deck	250 lbs.
	36"L x 80"W x 5"H	Flat Deck	350 lbs.

(1) Custom sizes available including 10" height and side bolsters

(2) Requires Accessory K-1365GFLX



MADE IN USA

1395 Pico Street, Corona, CA 92881
 Customer Service 1.866.KAP.MED.1 (888.527.6331)
 Facsimile 951.340.4361

Patented Product: US# 7454809, US# 8220090, US#7171711

Distributed By: