

# Bari™ 10A

## Bariatric Bed

**Medical**

Rental Solutions Delivered by



1 800 STRYKER Option 3

Universal | Hospital | Services®



**Overall aesthetic that allows the patient the dignity of a standard hospital bed look**

**BackSmart™ expansion with color-coded levers allow the bed to easily expand and retract to accommodate patients of all sizes**

**Centrally-located fifth wheel adds an additional pivot point, improving maneuverability and steering**

**Patient helper provides a BackSmart solution for patient repositioning**

**True 1000 lb. weight capacity**

### Large Capacity Patient Care

A practical solution for better patient care and comfort for all patients. Bari 10A can be used as a standard hospital bed, but can easily expand to accommodate your heaviest patients. Best-in-class expansion and retraction with color-coded levers expands 6 inches on each side. Centrally located fifth wheel improves maneuverability and steering.



Best-in-class expansion/retraction with color-coded levers.



Made with heavy-gauge, tubular steel, the Bariatric Patient Helper adds to patient dignity by allowing the patient to lift the upper body with arm strength.



A center-mounted fifth wheel provides an additional pivot point to turn the bed in a tighter arc. The fifth wheel retracts to clear the under-bed space.

# Bari 10A

## Bariatric Bed

**stryker**<sup>®</sup>



### Medical

3800 E. Centre Ave.  
Portage, MI 49002 USA  
t: 269 329 2100 f: 269 329 2311  
toll free: 800 787 9537

www.stryker.com

### Standard Features

- 1000 lb. patient capacity
- Fully articulated mattress support platform with independent back and knee movement, up/down adjustments, Trendelenburg and Reverse Trendelenburg positioning
- 36" wide mattress support platform sections easily expand to 48" wide
- 80" (expandable to 88") mattress support platform. Integrated foot end extension section extends an extra 8" length for taller residents/patients
- 18.75" to 33.75"
- Integrated scale system for resident/patient weight monitoring up to 1000 lbs
- Head and foot end siderails/assist devices
- One-handed siderail release
- Siderails on head and foot sections fold down for easy patient access
- Pendant control for all bed functions with separate selective lock-out to override functions
- Emergency manual release for mattress support platform
- Quiet DC motors
- Welded steel tubular frame
- Pan-style mattress support platform
- Wide stable base
- Heavy-duty casters
- Centerline caster locks
- Bed angle indicators
- Soft-tone frame color
- Head of bed indicator

### Specifications

• <b>Model Number:</b>	2230
• <b>Overall Length:</b>	Minimum: 90", Maximum: 98"
• <b>Overall Width:</b>	Minimum: 39", Maximum: 52"
• <b>Weight Capacity:</b>	1050 lbs.
• <b>Bed Frame Weight:</b>	700 lbs
• <b>Standard Height Range</b>	
High	33.7"
Low	18.7"
• <b>Litter Positioning</b>	
Maximum back angle to horizontal	65°
Maximum knee angle to horizontal	35°
Maximum knee to leg angle	45°
Trend	14.5°
Reverse Trend	15.3°
• <b>Patient Surface</b>	
Minimum	36" x 80"
Maximum	48" x 88"
• <b>Electronics</b>	
<b>Agency Approvals</b>	c-us UL, UL60601-1, IEC 60601-1-2, CAN/CSA C22.2 No. 601.1 IEC 60601-2-38, CAN/CSA C22.2 No. 601.2.38
<b>Volts</b>	120 Vac
<b>Ampere Rating</b>	5.8 amps
<b>Frequency</b>	60 Hz

*Stryker reserves the right to change specifications without notice.  
The Bari 10-A Bariatric Bed is manufactured for Stryker Medical.*

### Optional Features

- Foam surface
  - BariMatt 1000
  - BariMatt 1000 Plus
- Low air-loss surface
  - Arise 1000 EX
- IV rod
- Floor mat
- Patient Helper
- Fifth Wheel

Copyright © 2007 Stryker  
PC/SCG -7.5M-0508, 07-SMD-1938  
Mkt-Lit-198-092007 Rev D



**Mixed Sources**  
Product group from well-managed  
forests, controlled sources and  
recycled wood or fibre  
Cert no. SW-COC-1530  
www.fsc.org  
© 1996 Forest Stewardship Council