


## Quartet™ Mattress



- ▶ Three mattress models available with the following therapies:
  - Auto immersion
  - Low air-loss
  - Pulsation
  - Alternating pressure, turn assist or continuous lateral rotation
- ▶ Active sensors to monitor and automatically adjust immersion settings to match patient weight and position
- ▶ Internal dual battery for power outages and patient transfers
- ▶ Standard, bariatric, expandable and recessed deck mattress options

### The Quartet™ Mattress Series: Customizing the care your patients need

The OnCare® Quartet™ solves the challenge of managing high-acuity, complex wounds in one easy-to-use, comprehensive package. Engineered with intelligent auto immersion software that continuously monitors and adjusts to a patient's changing positions, the Quartet also offers low air-loss, pulsation and your choice of alternating pressure, turn assist or continuous lateral rotation. These features not only accelerate the recovery process, but reduce the burden on clinicians to continuously monitor patient movements and mattress settings, allowing for more time to spend on bedside care.





# Quartet™ Mattress Features and Specifications



## Clinical Applications:

- Stage I through IV pressure ulcers
- High-risk patients requiring aggressive wound prevention
- Excessive maceration
- Flaps and grafts
- Burns
- Acute surgical wounds

Therapy Features	Quartet Mattress Models		
	Auto Immersion	Turn Assist	Clear Turn
Auto Immersion	✓	✓	✓
Low Air-Loss	✓	✓	✓
Pulsation	✓	✓	✓
Alternating Pressure	✓		
Heel Zone Regulator	✓		
Turn Assist		✓	
Continuous Lateral Rotation			✓

## Standard Features for All Models

- Max Inflate
- Internal Dual Batteries
- Lock Out
- Alarms
- CPR Deflate
- Five Zone Mattress
- Fowler Boost
- Custom Comfort Settings
- Bottom Out Protection
- Breathable Non-Shear Nylon Top Sheet
- Cable Management System

Type	Quartet Mattress Specifications			
	Length	Width	Height	Weight Capacity
Standard	80"	36"	10"	Up to 500 lbs.
Bariatric	80" or 84"	48"	10"	Up to 1000 lbs.
Expandable	80"/84"/88"	36"/42"/48"	10"	Up to 1000 lbs.
Recessed	86"	36"	6" to 9"	Up to 500 lbs.

Visit [www.oncaremedical.com](http://www.oncaremedical.com) to learn more.

OnCare® Quartet™ is a proud offering of Agiliti Health, Inc. Patented Product US #8102270, US #7454809, US #8220090 and US #7171711.

### CORPORATE OFFICE

6625 West 78th Street, Suite 300, Minneapolis, MN, 55439

Toll Free: 800.847.7368 Phone: 952.893.3200 Fax: 952.893.0704 [www.agilitihealth.com](http://www.agilitihealth.com)

©2018 Agiliti Health, Inc. QUARTET:1218

