



To Prevent and Treat Pressure Ulcers

Plexus Aire Select Series

- Choice of mattress or overlay
- Mattress is effective for all stages of pressure ulcers
- Overlay is effective for prevention and early treatment of pressure ulcers
- Deep-cell, low-air-loss, selectable alternating pressure therapy
- User-friendly controller allows for customized therapy
- Optional foam safety bolsters
- Oversized top sheet cover for additional protection against shear available in nylon or Gaymar Silver, which is embedded with silver ions
- Suitable for all healthcare environments

Plexus® AIRE SELECT™

Low-Air-Loss/Alternating Pressure Mattresses and Overlays

Model Number P2500 and PX3580 Series

□ U.S.: 1.800.828.7341 716.662.2551 Canada: 1.800.203.4293 □ www.gaymar.com

 **GAYMAR®**
Pressure Ulcer Management

Plexus® Aire Select™

Dynamic Low-Air-Loss Therapy with Selectable Alternating Pressure Mattress and Overlays

The Control Unit:

- Effective low-air-loss to assist with moisture control
- Selectable Alternating Pressure with a five minute cycle time
- **Max. Inflate** for easy patient positioning
- Patient settings are restored after power outages
- CPR quick release provides for rapid deflation
- Control Unit is compact and lightweight

The Therapeutic Support Surface:

- Deep air cells constructed of high quality coated nylon for low shear and effective therapy
- Optional overlay (P200LAL25) provides pressure redistribution, alternating pressure, and low-air-loss in a convenient effective mattress overlay
- Plexus quilted top sheet is highly vapor permeable and oversized providing a low friction, low shear barrier. Machine washable and available in nylon or silver
- Unique modular design facilitates cleaning and interchangeability
- Convoluted foam base provides additional safety for your patient



P200LAL25 Overlay

Foot Stop/ZIP™ Cushion (Optional):

- Interchangeable Foot Stop/ZIP cushion permits tailoring the mattress length to fit the patient while reducing pressure in the critical heel area
- ZIP cushion can be independently used anywhere along the mattress

Ordering Information:

Catalog #	Description	Dimensions	Weight
C2500	Control Unit	10" w x 12" h x 5" d (25cm x 30cm x 13cm)	9.5lb (4.3kg)
M2500	Aire Select Mattress	34" w x 80" l x 10" h (86cm x 203cm x 25cm) <i>Also available in widths of 32", 39", 42", 48", 54", and 60" (81, 99, 107, 121, 137, 152cm)</i>	15lb (6.8kg)
M200LAL25	Aire Select Overlay	34" w x 76" l x 3.5-6" d (86cm x 193cm x 9-15cm)	10lb (4.5kg)
X3580	Aire Select Safety Mattress	35" w x 80" l x 10" h (89cm x 203cm x 25cm) with 13" (33cm) side bolsters	27lb (12.2kg)
P2500	Aire Select Mattress Replacement System: Mattress with Control Unit		
P200LAL25	Aire Select Overlay System: Overlay with Control Unit		
PX3580	Safety Mattress System: Safety Mattress with Control Unit		

Specifications:

Electrical Rating:	120V~, 50Hz, 5A
Agency Approvals:	UL 544, EN 60601-1-2
Warranty:	Control Unit: 2 year limited, Mattress/Overlay: 1 year limited
Weight Limit:	350 lbs.
Medicare HCPCS Code:	E0277

Aire Select Safety Mattress

- 13" (33cm) firm foam side bolsters with or without egress section
- Zipped cover with quilted topper available in nylon or silver, which is embedded with silver ions



Optional Foot Stop/ZIP cushion uniquely maintains proper patient positioning while *reducing* interface pressure on the critical heel area.



Separate carrying case also available.



Gaymar Industries, Inc.

10 Centre Drive

Orchard Park, NY 14127

Toll Free 1.800.828.7341

Fax 1.800.993.7890

International +1 716.662.8636

International Fax +1 716.662.0730

www.gaymar.com

Available in Canada through

Waterloo Bedding Co. Ltd.

825 Trillium Drive

Kitchener, ON N2R 1J9

Toll Free: 1.800.203.4293 • 519.742.4447

Fax: 519.742.3699