

- seat configurations from 20" wide to 32" wide and 18" deep to 22" deep for client comfort, positioning and custom fit
- designed for clients weighing from 700 to 850 pounds
- folding frame construction for portability
- exclusive double square cross braces for rugged durability
- adult, hemi and super hemi seat to floor height
- removable or optional height adjustable removable arms
- height adjustable back for a custom fit
- all seat configurations available with power assist reclining back

Regency XL 2000



Bariatric Rehab Wheelchairs

Model	Seat W x D	Back Style	Arm Style	Seat to Floor Height*	Reclining Back	Weight Cap.**
6720-18	20 X 18	16" 17" 18" height adj.	removable	15" 16" 17" 18"	optional	700 to 750
6722-18	22 X 18	16" 17" 18" height adj.	removable	15" 16" 17" 18"	optional	700 to 750
6724-18	24 X 18	16" 17" 18" height adj.	removable	15" 16" 17" 18"	optional	700 to 750
6726-18	26 X 18	16" 17" 18" height adj.	removable	15" 16" 17" 18"	optional	700 to 750
6728-18	28 X 18	16" 17" 18" height adj.	removable	15" 16" 17" 18"	optional	700 to 750
6730-18	30 X 18	16" 17" 18" height adj.	removable	15" 16" 17" 18"	optional	700 to 750
6732-18	32 X 18	16" 17" 18" height adj.	removable	15" 16" 17" 18"	optional	700 to 750
6720-20	20 X 20	16" 17" 18" height adj.	removable	15" 16" 17" 18"	optional	700 to 750
6722-20	22 X 20	16" 17" 18" height adj.	removable	15" 16" 17" 18"	optional	700 to 750
6724-20	24 X 20	16" 17" 18" height adj.	removable	15" 16" 17" 18"	optional	700 to 750
6726-20	26 X 20	16" 17" 18" height adj.	removable	15" 16" 17" 18"	optional	700 to 750
6728-20	28 X 20	16" 17" 18" height adj.	removable	15" 16" 17" 18"	optional	700 to 750
6730-20	30 X 20	16" 17" 18" height adj.	removable	15" 16" 17" 18"	optional	700 to 750
6732-20	32 X 20	16" 17" 18" height adj.	removable	15" 16" 17" 18"	optional	700 to 750
6720-22	20 X 22	16" 17" 18" height adj.	removable	15" 16" 17" 18"	optional	700 to 750
6722-22	22 X 22	16" 17" 18" height adj.	removable	15" 16" 17" 18"	optional	700 to 750
6724-22	24 X 22	16" 17" 18" height adj.	removable	15" 16" 17" 18"	optional	700 to 750
6726-22	26 X 22	16" 17" 18" height adj.	removable	15" 16" 17" 18"	optional	700 to 750
6728-22	28 X 22	16" 17" 18" height adj.	removable	15" 16" 17" 18"	optional	700 to 750
6730-22	30 X 22	16" 17" 18" height adj.	removable	15" 16" 17" 18"	optional	700 to 750
6732-22	32 X 22	16" 17" 18" height adj.	removable	15" 16" 17" 18"	optional	700 to 750
6826-20	26 X 20	16" 17" 18" height adj.	removable	15" 16" 17" 18"	optional	850
6828-20	28 X 20	16" 17" 18" height adj.	removable	15" 16" 17" 18"	optional	850
6830-20	30 X 20	16" 17" 18" height adj.	removable	15" 16" 17" 18"	optional	850
6832-20	32 X 20	16" 17" 18" height adj.	removable	15" 16" 17" 18"	optional	850
6826-22	26 X 22	16" 17" 18" height adj.	removable	15" 16" 17" 18"	optional	850
6828-22	28 X 22	16" 17" 18" height adj.	removable	15" 16" 17" 18"	optional	850
6830-22	30 X 22	16" 17" 18" height adj.	removable	15" 16" 17" 18"	optional	850
6832-22	32 X 22	16" 17" 18" height adj.	removable	15" 16" 17" 18"	optional	850

Note*: 15" seat to floor height requires 20" rear wheels and 6" caster wheels

Note**: 750 pound capacity requires wheel package up-grade, see price list for details

For more information or to arrange for a free demonstration of the Bariatric Rehab Wheelchairs please call toll free


1-800-537-2521

GENDRON 

400 East Lugbill Rd.
Archbold, OH 43502
Fax 419-446-2631
www.gendroninc.com
email sales@gendroninc.com



Removable and Height Adjustable Arm Rest



Rear Wheel Positions for Seat Heights From 15" to 18"

3-Position Adjustable Back Height